


Article Title	A Love Hate Relationship with Flowable Composites
Category	<input type="checkbox"/> Dental Adhesives <input checked="" type="checkbox"/> Composites <input type="checkbox"/> Cements <input type="checkbox"/> Impressioning materials <input type="checkbox"/> Prevention
Keyword(s)	Filtek composite, flowable composite
Hero Visual (in high resolution also as separate file)	
Copyright image	3M
Alt text image	Love for flowable composites
Meta description (<150 characters)	We love flowable composites because of the ease of placement and the excellent wetting of the tooth with its lower viscosity. We are tormented, by the constant threat of bubbles.
Author	Dr. David Clark, DDS About Dr. David Clark 25 years of clinical microscope dentistry in a hometown practice, lab research, and hundreds of lectures and hands on courses has enabled him to produce matrix systems and methods that are reinventing, not just the way we do restorative dentistry, but more importantly, the way that we think about restorative dentistry.
Text (>750 words) in web format: <ul style="list-style-type: none"> • Internal links • Precise Product Names & urls from the 3 M 	<p>We love flowable composites because of the ease of placement and the excellent wetting of the tooth with its lower viscosity. We are tormented, however, by the constant threat of bubbles that so often appear as we syringe the material.</p> <p>Let's quickly review the current role of flowable composites.</p> <p>When the first dental resin-based composites were first</p>

- product catalogue
- Bullet points style
- Imagery
-

Infographics

introduced by 3M in 1964, they were thought of as an “amalgam substitute”. In other words, we thought that we needed to place it in little increments and pack it, sculpt it, and carve it like we did with amalgam.

When the first flowable composites were introduced in 1996,¹ nobody was quite sure what to do with them, but we could see the definite advantages of an injectable material that instantly wets the tooth structure. What we have learned over time is that the more we manipulate composites, the more problems we can experience. Composite today is moving to being injected with little or no hand manipulation and flowable is an ideal injectable material.

The trouble with bubbles ...

The biggest complaint that clinicians report to me regarding flowable composites, is the bubbles that appear during syringing of the flowables. These **bubbles in the flowable** can be maddeningly difficult to eliminate. The bubbles are difficult to “pop” because of the viscosity of the resin. In desperation we give up on popping them and try to drag them outside of the cavity prep.



Two different flowables have been syringed onto dentoform teeth. Both samples have multiple bubbles. The bubbles in the lighter shade (right) are more difficult to see.

Two times when bubbles are discovered too late ...

1. Bubbles visible on radiographs



In follow-up radiographs, two bubbles are visible for this endodontically treated tooth. It is embarrassing and hard to explain to patients. These bubbles drive endodontists crazy.

2. Blemishes on the restoration caused by bubbles

When a bubble leaves a divot on the surface of a composite, blemishes with embarrassing discoloration can occur.



A small bubble in the composite has significant stain that the patient complained about, and in fact was the reason they left their previous dentist.



This bubble in the conservative Class I Fissurotomy preparation

	<p>is an all too common sight.</p> <p>Conclusion</p> <p>There's a lot to love about flowable composites. They've become integral to the modern practice for most restorative dentists. While flowable materials and modern techniques continue to improve, I also see promise in new delivery systems. I've been evaluating 3M's new syringe design for their flowable composites. This new innovative design has the potential to reduce bubbles in our procedures and give dentists more control.</p> <p>Related 3M products:</p> <p>For aesthetics, choose 3M™ Filtek™ Supreme Flowable Restorative (ADD PRODUCT DETAIL PAGE URL)</p> <ul style="list-style-type: none"> • Excellent adaptation, polish retention and wear resistance • Features 3M's TRUE nanotechnology • Class III and V restorations up to 2mm • Base/liner under direct restorations <p>For efficiency, use 3M™ Filtek™ Bulk Fill Flowable Restorative (ADD PRODUCT DETAIL PAGE URL)</p> <ul style="list-style-type: none"> • Excellent adaptation, low shrinkage • Ideal as a liner, bulk fill base, sealant and more • Class III and V restorations up to 4mm • Base/liner under direct restorations
Recommended products	<ul style="list-style-type: none"> • 3M™ Filtek™ Supreme Flowable Restorative • 3M™ Filtek™ Bulk Fill Flowable Restorative
Sources Footnotes	<p>1. 1 https://pdfs.semanticscholar.org/8de5/add978297f750ce7dc83aeb3c3ce1a2648a0.pdf</p> <p>Original article: Our Love/Hate Relationship with Flowable Composites - 3M Dental Blog</p>
Supportive imagery, Infographics, etc.	
Related articles	